


# Information



The following information merely states the requirements for the computer system. Even if the stated system requirements are met, correct operation of the DÜRR MEDICAL hardware/software can still be disrupted by specific hardware and software features on the customer side. In these cases, DÜRR MEDICAL assumes no liability for ensuring the trouble-free operation of the DÜRR MEDICAL hardware/software.

The connection of additional systems to the PC may change the system requirements. Comply with the system requirements of all connected systems.



## Image plate scanner

The system requirements of the software being used to operate the device must always be met during its operation.

We recommend using the latest version of the software.

If third-party software is being used to operate the device, compliance with its system requirements must be assured.

The following additional system requirements must also be met:

### CR 7 VET

Interface:	USB 2.0 USB 3.0 Ethernet $\geq 100$ MBit/s
Total maximum cable:	USB: up to 5 m with USB extension cable USB: up to 20 m with additional USB repeater / active USB hub (max. 5 m each) - details available on request
Software:	Vet-Exam Pro version 1.0 or higher Vet-Exam Pro Connect version 1.2 or higher Vet-Exam plus version 7 or higher

### CR 8 VET

Interface:	USB 2.0 USB 3.0 Ethernet $\geq 100$ MBit/s WLAN IEEE 802.11b/g with WPA and WPA2 encryption
Graphics card:	Resolution $\geq 1024 \times 768$ Depth of colour 32-bit, 16.7 million colours
Software:	Vet-Exam Pro version 1.0 or higher Vet-Exam Pro Connect version 1.2 or higher

### CR 35 VETwin

Interface:	Ethernet $\geq 100$ MBit/s WLAN IEEE 802.11b/g with WPA and WPA2 encryption
Software:	Vet-Exam Pro version 1.0 or higher Vet-Exam Pro Connect version 1.2 or higher Vet-Exam plus version 7 or higher



**DÜRR NDT GmbH & Co. KG**

Division DÜRR MEDICAL

Höpfigheimer Str. 22

74321 Bietigheim-Bissingen, Germany

[www.duerr-medical.de](http://www.duerr-medical.de)

**9000-608-100/02**



**DÜRR  
MEDICAL**

2204V003

Small X-ray system

The system requirements of the software being used to operate the device must always be met during its operation. We recommend using the latest version of the software. If third-party software is being used to operate the device, compliance with its system requirements must be assured. The following additional system requirements must also be met:

VistaRay 7

Interface:	USB 2.0 USB 3.0
Total maximum cable:	USB: 5 m
Software:	Vet-Exam Pro version 1.0 or higher Vet-Exam Pro Connect version 1.2 or higher Vet-Exam plus version 7 or higher

Intraoral camera systems

The system requirements of the software being used to operate the device must always be met during its operation. We recommend using the latest version of the software. If third-party software is being used to operate the device, compliance with its system requirements must be assured. The following additional system requirements must also be met:

VistaCam iX HD Smart

Interface:	USB 2.0 USB 3.0
Total maximum cable length:	Up to 5 m with USB extension cable Up to 20 m with additional USB repeater / active USB hub (max. 5 m each) - details available on request
Software:	Vet-Exam Pro version 1.0 or higher Vet-Exam Pro Connect version 1.2 or higher

VistaCam iX

Interface:	USB 2.0 USB 3.0
Total maximum cable length:	Up to 5 m with USB extension cable Up to 20 m with additional USB repeater / active USB hub (max. 5 m each) - details available on request
Software:	Vet-Exam Pro version 1.0 or higher Vet-Exam Pro Connect version 1.2 or higher

## Software

### Vet-Exam Pro

Valid from version 1.2

CPU:	≥ Intel Core i3
RAM:	≥ 4 GB
Operating systems:	64-bit operating system: Microsoft Windows 10 (Pro or higher) Microsoft Windows 11 (Pro or higher) Microsoft Windows Server 2019 Microsoft Windows Server 2022
Hard disk:	Workstation (without database) ≥ 50 GB The database memory requirements depend on the number of images taken at the surgery in question. (Camera image: approx. 1 MB / X-ray image: approx. 2 MB - 10 MB)
Data backup:	Daily data back-up
Interface:	Ethernet ≥ 100 MBit/s
Graphics card:	Resolution ≥ 1280 x 1024 Depth of colour 32-bit, 16.7 million colours

### Vet-Exam plus

Valid from version 7

CPU:	≥ Intel Pentium IV compatible, 1.4 GHz
RAM:	≥ 1 GB (2 GB recommended))
Operating systems:	Microsoft Windows 10 (Pro or higher) Microsoft Windows Server 2016 Microsoft Windows Server 2019
Hard disk:	Workstation (without database) ≥ 50 GB The database memory requirements depend on the number of images taken at the surgery in question. (Camera image: approx. 1 MB / X-ray image: approx. 2 MB - 10 MB)
Drive:	DVD-ROM
Data backup:	Daily data back-up
Interface:	Ethernet ≥ 100 MBit/s
Graphics card:	Resolution ≥ 1024 X 768 Depth of colour 32-bit, 16.7 million colours